

## KN1540BA – Aquasep



Anti-bacterial solution for portable eye washes stations with a volume of 4 gallons to 20 gallons (15litres to 75litres). Including Krusman models KN1540 (60litres) and KN1534 (34litres).

8 Fluid ounces (236mls) should be added to a minimum water level volume of 4 gallons (15litres) and up to a maximum of 20 gallons (75litres). Contact sales representative for MSDS document.

It is recommended that the station be cleaned and the water replaced at a minimum of every 6 months. Examination at weekly intervals which give evidence of any particulates or discolouration suggests the water be replaced immediately.

If you have any further questions please contact your local sales representative or email: [info@krusman.se](mailto:info@krusman.se)